

NE[®]



Code: NE026 a/b/c

SPONSTIM

tCS-only

The stimulation
sponge electrode



NEUROELECTRICS
BARCELONA
AV. Tibidabo 47 bis
08035 Barcelona
Spain



DESCRIPTION

The Sponstim consists of a sponge cover, a carbon rubber core and a metallic pin made of nickel plated brass.

The contact surface differs among the three models:

- (a) circular shape with 25 cm²
- (b) circular shape with 8 cm²
- (c) rectangular shape with 35 cm²

MAINTENANCE

After each use, wash the sponges with tap water and let them dry before storage. When stored, make sure the metallic pins do not come into contact with the sponges, so they do not become rusty.

DURABILITY

The cleaner Sponstim is kept, the longer it lasts. On average, after 100 hours of stimulation, the rubber core loses its conductive properties, and the electrode should be replaced.

