


NE[®]



Code: NE027

EARCLIP

Reference

The dual reference
CMS & DRL electrode



NEUROELECTRICS
BARCELONA
AV. Tibidabo 47 bis
08035 Barcelona
Spain



The electrical reference Earclip is a dual electrode used to connect both CMS and DRL channels simultaneously to the same earlobe. Compared to the Sticktrode, the Earclip is a re-usable alternative that offers a quick assembling process.

DESCRIPTION

The Earclip consists of two opposed Ag/AgCl pellets of 8 mm diameter in a clip. Each one of the two pellets has a 0.5 cm² contact surface.

MAINTENANCE

After each use, clean the Earclip using a tissue. If needed, you may use water to complete the cleaning. Avoid direct sunlight exposure and contact with metals when the electrode is not being used.

DURABILITY

The Earclip is a long lasting electrode so it can be reused many times. Its intended number of uses is 100.

USE INSTRUCTIONS

1 – Apply some electrode gel to the sensors (recommended).

2 – Connect the CMS and DRL cable connections to the Earclip.

3 – Clip it to the right or left earlobe.

